

Experience a new  
way of handling  
your mobile lift.

**EVADrive**  
System **ROMedic**



handicare

# Why we created EvaDrive.

In healthcare the cost of manual handling incidents, especially back injuries is constantly increasing as is the number of people requiring care and transfer assistance. Further, inadequate budgets in combination with the onset of reduced care packages and an increasing population of "plus size" patients add to the on going risk and potential for load injuries.

Our goal is to provide new and innovative lifting and transfer aids that reduce the risk to carers and avoid these load injuries. That is why we created EvaDrive.

EvaDrive is a new motorized mobile lift which greatly reduces the force required to maneuver the patient while in the lift. EvaDrive runs almost effortless, even in narrow spaces. With its intuitive control handle it is as easy to operate as your ordinary shopping trolley!

When used with Handicare's functional and easy-to-use SystemRoMedic™ slings and accessories, the EvaDrive is the right choice for lifting in connection with transfers to and from the bed, chair or floor as well as for horizontal repositioning and gait training.





*Intuitive control panel*

*Perceptive handle*

*Powerful lithium battery*

*Turnable rear wheels*

# Intuitive operation. A completely new way to maneuver your mobile lift.

EvaDrive operates intuitively and almost without force. The handle perceives your every move. With the tip of your finger you can move EvaDrive into very narrow spaces, when for example repositioning a patient in bed in a small hospital room.

EvaDrive will turn, reverse, and run forward simply by gently touching the handle in the desired direction. EvaDrive knows where you want to go – and takes you there!



A single carer can easily operate the lift with a plus size user without using any force.





handicare

EVA 450FE

EVA DRIVE

# How it works.

## The control panel

---



*Press the On/Off button and EvaDrive is ready to operate within seconds.*



*The rear wheels can be turned to any angle.*



*When you need high precision and low speed - simply press the snail.*



*The service indicator will announce when EvaDrive is due for inspection.*

## The handle

---



*The perceptive and ergonomic handle is equipped with load cells that intercepts your every move. By gently pushing it, EvaDrive moves forward. By pulling it lightly, EvaDrive reverses. The handle will sense your every desired movement. The EvaDrive will go in any direction at the point of your fingertips. The handle is ergonomic and adjustable in height. Also, a safety switch will prevent any undesired movement.*

## The wheels

---



*The rear wheels turn to any angle by using the control panel buttons. EvaDrive moves into very narrow spaces enhanced by the lateral movement of the wheels.*



*EvaDrive easily operates in narrow spaces.*

## The battery

---

*The powerful lithium battery will last an entire work shift. Charge the battery when EvaDrive is not in use, in order to have a fully operational lift ready for the next work shift. The battery is detachable and can be replaced if needed.*



