




WendyLett4Way with handles

SystemRoMedic™

Trial fitting protocol

Patient's name: _____

Social insurance number: _____



Purpose of assistive device:

When manual positioning in bed requires too much physical exertion, when there are not enough caregivers, or when the patient is too heavy or is passive, WendyLett4Way with handles can be used in combination with a lift. This combination helps to reduce friction under the patient and the caregiver can care for the patient more easily in bed.

Patient group:

Patients with e.g., pain, anxiety, stiffness after surgery, obesity, spasticity or low muscle tone.

In which situations is the assistive device used?

- | | |
|---|--|
| <input type="checkbox"/> Occasional turning to right side | <input type="checkbox"/> Occasional turning to left side |
| <input type="checkbox"/> Positioning on right side | <input type="checkbox"/> Positioning on left side |
| <input type="checkbox"/> Hygiene care | <input type="checkbox"/> Sore dressing |
| <input type="checkbox"/> Higher up in bed | <input type="checkbox"/> Dressing and undressing |
| <input type="checkbox"/> Application of sling | |
| <input type="checkbox"/> Misc: _____ | |

Assistive device used with WendyLett4Way with handles (1657):

- WendyLett slide sheet 140x200 cm (1634) WendyLett slide sheet 90x200 cm (1641)

Complete package:

- WendyLett complete package 1628 (1634x1657)
 WendyLett complete package 1646 (1641x1657)
 WendyLett complete package 1656 (1651x1657)

Accessories that can be used with this product:

- AdjustmentLoop (extension strap) - 18 cm (49182006)



Lyft

- Two-point slingbar, width 45 cm
 O verhead lift with two-point slingbar, traverse (recommended) / straight rail
 Mobile lift with med two-point slingbar

Preparations/Positioning in bed:

Instructions for use

Patient's name: _____

Social insurance number: _____

Which size of WendyLett is used?

WendyLett 140x200 cm

WendyLett 200x200 cm alt. 200x90 cm with lace lock

Has the caregiver using this assistive device:

Read the instruction guide?

Yes

No

Received practical instruction in the transfer of patients?

Yes

No

Instructions for caregivers:

E.g., Raise the gate on the side you are facing. Positioning of head/arms/hands? Positioning of legs/feet? Are bed functions used?

Number of caregivers: _____

Instructions for patient:

E.g., In what way can the patient be active/participate? Which instructions are given?

Trial-fitted by occupational therapist/physiotherapist:

Date: _____

Telephone: _____

Telephone: _____

Handicare AB
Veddestav. 15, Box 640
SE-175 27 Järfälla
SWEDEN

Tel: +46 (0)8-557 62 200
Fax: +46 (0)8-557 62 299
E-mail: info@handicare.se

Internet: www.handicare.com

Changes or questions, please contact the physiotherapist!

Manuals and instructions about function tests is available at the Handicare website
www.handicare.com

